

Ampera™ bed – Classic

An affordable care option



Features and benefits

The **lower height of the bed** assists in minimising the risk of fall-related injuries

A **220 kg Safe Working Load** accommodates a greater range of patients/residents

Central braking and steering functions ensure enhanced security for patients/residents

Line-Of-Site™ backrest indicator ensures accurate patient/resident positioning information for caregivers

Intelligent AutoContour™ and Sliding Backrest minimise patient/resident migration in bed, reducing abdominal compression during bed articulation

Caregiver control with total lock-out means enhanced operational safety

Ampera bed - Classic



Pendant control



Bilateral mechanical CPR



Central braking



Optional: Full-length wooden siderails

Technical Specifications (for standard configuration)

Safe Working Load	220 kg
Bed weight	137 kg
Height adjustment (dependent on castors)	
standard	35 to 73 cm
optional	39 to 77 cm
Protection class	IPX4
Bed angles	
Head section	0 to 65°
Thigh section	0 to 28°
Foot section	-3 to -22°
Sleep deck width	90 cm
Sleep deck length	200 cm
Sleep deck surface	metal
Siderails	
standard	none
optional	3/4 length metal
optional	full-length wooden

Castors

standard	4 x 100 mm single band
optional	4 x 125 mm single band
optional	4 x 150 mm single band

Steering castor

standard	at head end
optional	at foot end

Central brake

at foot end

Lifting/IV pole holder

inside double holder
at head end

Bilateral mechanical CPR (Cardiopulmonary Resuscitator)

Controls

see controls sheet

Designs/Colours

see designs and colours sheet

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

The Ampera medical bed is intended for use in residential care facilities and retirement homes for patients with a weight ≥ 40 kg, height ≥ 146 cm and BMI ≥ 17 .
Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 Zi du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN138302-04, 10 July 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom