

# Ampera™ bed – Style

For maximum versatility



## Features and benefits

The **lower height of the bed** assists in minimising the risk of fall-related injuries

A **220 kg Safe Working Load** accommodates a greater range of patients/residents

**Central braking and steering functions** ensure enhanced security for patients/residents

**Line-Of-Site™** backrest indicator ensures accurate patient/resident positioning information for caregivers

**Intelligent AutoContour™ and Sliding Backrest** minimise patient/resident migration in bed, reducing abdominal compression during bed articulation

**Caregiver control with selective bilateral lock-outs** means enhanced operational safety

**Patented bed height memorisation and one-touch egress positioning** allows patients/residents to select a preferred height and position for easier and safer bed exit

**Backup battery** provides continuous functioning capabilities

**Egress positioning** contributes to easier patient/resident mobility

**Automatic horizontal stop position** notifies caregivers when the patient/resident is in a flat position

## Ampera bed - Style



Pendant control



Bilateral mechanical CPR with selective lock-outs



Optional: Flexible arm control



Optional: Full-length wooden siderails

### Technical Specifications (for standard configuration)

Safe Working Load	220 kg
Bed weight	137 kg
Height adjustment (dependent on castors)	
standard	35 to 73 cm
optional	39 to 77 cm
Protection class	IPX4
Bed angles	
Head section	0 to 65°
Thigh section	0 to 28°
Foot section	-3 to -22°
Trendelenburg/Reverse Trendelenburg	-13 to +13°
Sleep deck width	90 cm
Sleep deck length	200 cm
Sleep deck surface	
standard	metal
optional	High Pressure Laminate removable (HPL)
Siderails	
standard	3/4 length metal
optional	full-length wooden
Siderail height	38 cm

### Castors

standard	4 x 100 mm single band
optional	4 x 125 mm single band
optional	4 x 150 mm single band

### Steering castor

standard	at head end
optional	at foot end

### Central brake

at foot end

### Lifting/IV pole holder

inside double holder  
at head end

### Bilateral Mechanical CPR (Cardiopulmonary Resuscitator)

### Service Required Indicator

### Controls

see controls sheet

### Designs/Colours

see designs and colours sheet

### Additional options

Egress handles at head section\*

(available only in combination with optional HPL sleep deck)

EasyChair™ positioning

\* Not available with full-length wooden siderails

Hill-Rom is a leading global medical technology company with more than 10,000 employees worldwide. We partner with healthcare providers in more than 100 countries, across all care settings, by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency and Respiratory Health.** The people, products and programs of Hill-Rom work towards one mission: **Every day, around the world, we enhance outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

[www.hill-rom.com](http://www.hill-rom.com)

The Ampera medical bed is intended for use in residential care facilities and retirement homes for patients with a weight  $\geq 40$  kg, height  $\geq 146$  cm and BMI  $\geq 17$ .

Class I

Manufacturer: Hill-Rom S.A.S. – BP 14 Zi du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.  
Doc. No: 5EN138303-05, 20 June 2017

Enhancing outcomes for patients and their caregivers:

**Hill-Rom**