



Hill-Rom® Procedural Stretcher

Go Beyond Mobility™

Customizable to be as simple or as versatile as you need, our Procedural Stretcher accommodates all patient conditions as well as your needs in this uncertain environment.

A rush to the operating room is not just about running down hallways and turning corners. It's about the minute details that contribute to a safe and comfortable experience for both you and your patient.

Our stretchers **provide more than mobility** for you and your patient:

- A safe and comfortable patient experience
- Simple and safe maneuvering
- Durability to withstand your demands

A Safe and Comfortable Patient Experience

Patients are spending more time on stretchers, and as a result, the patient experience, including pressure injury prevention, has become increasingly important.

- **Help prevent pressure injuries** through your choice of pressure redistribution surfaces.
- **Protect your patient's hands** with bumper protected siderails.
- **Prevent your patient from migrating** down our stretcher with Auto Contour™ feature.
- **Provide a gentle patient ride** with spring-loaded Steering Plus™ system.
- **Allow for safe patient ingress and egress** with low height.
- **Provide more coverage** for smaller patients with enhanced safety siderails.



Enhanced safety siderails



Simple and Safe Maneuvering

Whether you're making crisp turns down narrow halls, getting the stretcher rolling from a stop, or nestling it into a crowded elevator, our stretchers are ideal in any setting. Hill-Rom Stretchers are easy to maneuver to help reduce your risk of injuries.

- **Simple activation** requires less force to activate steer pedal than Stryker's Big Wheel.
- **360° Maneuvering** with Steering Plus™ Systems that grip the floor and reduce your turning radius.
- **Ease the burden on your back;** BackSaver Fowler® feature uses patient's weight to help position the head section.
- **Zero Transfer Gap** — OneStep® Tuck-Away Siderails create a foundation for seamless patient transfers.
- **Obtain an immediate weight without having to move your patient** with an integrated dual-display scale option that memorizes last weight.

Durability to Withstand your Demands

At Hill-Rom, we understand that because caregivers have to do more with less, having a reliable and durable product is more important than ever.

- **700 lb weight capacity** for all stretchers
- **Compatibility** with lifts and overbed tables
- **Easy to clean** smooth deck and siderail design
- **Power washable**

Go Beyond Mobility™

A safe and comfortable experience for both you and your patient starts on a Hill-Rom® stretcher.

A SAFE & COMFORTABLE PATIENT EXPERIENCE	STANDARD FEATURES	OPTIONAL FEATURES
26" width	X	
30" width		X
Auto Contour™ Mechanism (includes BackSaver® and Knee Gatch)		X
Manual Knee Gatch		X
Bumper Protected Siderails	X	
Enhanced Safety Siderail		X
Trendelenburg/Reverse Trendelenburg	X	
Six Transport Strap Locations	X	
PRESSURE REDISTRIBUTION MATTRESS OPTIONS		
3" Comfortline® foam mattress	X	
4" or 5" Comfortline® foam mattress		X
5" AccuMax Quantum™ mattress (Elastimax or Hypolex)		X
SIMPLE & SAFE MANEUVERING		
Steering Plus™ System	X	
OneStep® Tuck-Away Siderails	X	
BackSaver® Fowler feature		X
Brake/Steer – head and foot	X	
Brake/Steer – all four sides		X
Dual-Sided hi/lo Hydraulic Control Pedals	X	
3-Sided hi/lo Hydraulic Control Pedals		X
Ergonomic Push Handles		X
Ergonomic Push Handles with Integrated IV Pole Transport Holder		X
Active Hand Brake		X
Integrated O ₂ tank holder with utility tray		X
3/4" Standard IV Pole Sockets	X	
Removable IV Pole – 2 Stage		X
Permanent IV Pole – 3 Stage (40 lb weight capacity)		X
Integrated Dual-Display Scale		X
Radiolucent Deck		X
Upright Chest X-Ray Cassette Holder		X
DURABILITY TO WITHSTAND YOUR DEMANDS		
700 lb weight capacity	X	
Full Perimeter Bumpeting	X	
Smooth Deck Design	X	
Decals Options	X	
Bumper Color Options – choice of neutral, blue, red, teal or purple (green bumpers only for IntelliDrive® models)	X	

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
171254 rev 6 29-AUG-2018 ENG – US

For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com

Technical Specifications

Maximum weight limit	700 lbs (317.5 kg)
Low position (floor to top of sleep deck)*	20.7" (52.6 cm)
High position (floor to top of sleep deck)	34.25" (87 cm)
Overall length	83" (210.8 cm)
Overall maximum width (siderails up)	31" (78.7 cm) or 35.5" (90.2 cm)
Overall maximum width (siderails stored)	29.25" (74.3 cm) or 33.25" (84.5 cm)
Siderail length	47" (119.4 cm)
Siderail height above sleep deck	14" (35.6 cm)
Comfortline® mattress size	26" (66 cm) x 75" (190.5 cm) or 30" (76.2 cm) x 75" (190.5 cm)
AccuMax® VPC mattress size	26" (66 cm) x 76" (193 cm) or 30" (66 cm) x 76" (193 cm)
Maximum head elevation	90° 80° with BackSaver Fowler™ or Auto Contour™
Maximum Trend/Reverse Trend	18°
Floor to base clearance	3.5" (8.9 cm)
Below IntelliDrive® Power Drive module	1.5" (3.8 cm)
Caster size	8" (20.3 cm)
Electrical	120V AC, 60 Hz

Warranty

Frame & Welds for life of product
Three-year parts
Two-year mattress
One-year labor

Accessory Options — Please confirm accessory availability for model of stretcher ordered

Accessory IV pole	Liquid oxygen tank holder
Ankle stirrups	Oxygen tank holder
Armboard	Paper roll dispenser
Chart holder	Patient tray
Convertible footboard	Restraint straps
Detachable footboard	Siderail covers
Foot extender	Stretcher pillow
Infusion support system	Transport straps
IV transporter	Utility tray
Lateral X-ray cassette holder	

*Dependent on configuration

Enhancing outcomes for patients and their caregivers:

Hill-Rom