




Hillrom™



TS 7500 Surgical Table Platform with floor-mounted column and IMRIS' MR Conditional MR Neuro Tabletop.


MR NEURO OR TABLETOP FOR IMRIS® SURGICAL THEATRE

SPECIALTY OR TABLETOP FOR NEUROSURGERY WITH TS 7500 OR TABLE

Enhanced Workflow for IMRIS Surgical Theatre

Support complex neurosurgery with Hillrom's MR Neuro OR Tabletop.

- Optimal surgical site access due to individual two-segmented tabletop.
- Anatomical adjusted head positioning with special head fixation devices.
- Ideal intraoperative imaging access due to lean design.
- Enhanced flexibility with interchangeable tabletops and the TS 7500 Surgical Table Platform.



**APPLICATION:
NEUROSURGERY**



PRECISE & FLEXIBLE



OPTIMAL INTRAOPERATIVE ACCESS

Adaptable Head Fixation Devices (HFD) 100/200 by IMRIS for flexible head positioning.



LIMITLESS FLEXIBILITY

TS 7500 MR Neuro OR Table and IMRIS® iMRI System offer various positioning options.



CUSTOMIZED WORKFLOW

Enhance patient safety and increase workflow with scanning position memorization.



GENERAL TO SPECIALTY

Transform your OR to accommodate multiple specialties with one single column for the TS 7500 OR Table and a wide range of tabletops.

TECHNICAL DETAILS

Material no.	# 2064888
Dimensions (L x W)	95" x 20.3" - 21.3" (2412 mm x 515 mm 540 mm)
Height Range	29.8" - 51.5" (758 mm - 1308 mm)
X-ray area	Head Support
Longitudinal slide	18.6" (472 mm)
Trendelenburg / reverse Trendelenburg	± 45°
Tilt	± 30°
Maximum patient weight	397 lbs (180 kg)

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

IMRIS is a registered trademark of DEERFIELD IMAGING, S.À R.L

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.