



TS 7500 OR Table with floor mounted column and U24 tabletop;
TS 7500 OR Table wireless remote (below).

UNIVERSAL TABLETOPS

SPECIALTY SURGICAL TABLETOPS FOR A WIDE RANGE OF PROCEDURES

The TS 7500 Operating Table offers a high number of positioning options for a wide variety of applications: from general and neurosurgery to gynecology, urology, traumatology, and orthopedics. Flexibility is further increased by an extensive range of positioning pads and special table accessories.



**BARIATRIC****GENERAL SURGERY****NEUROSURGERY****SPINAL SURGERY****THORACIC SURGERY****ORTHOPEDICS****GYNECOLOGY****ORTHOPEDICS AND TRAUMA****TRANSFER TABLETOPS**

TECHNICAL DETAILS

	U14 Tabletop	U24 Tabletop	U26 Tabletop
Material no.	# 1846376	# 1846381	# 1846386
Dimensions (L x W)	32.8" x 23.6" (832 mm x 600 mm)	42.1" x 23.6" (1070 mm x 600 mm)	42.1" x 23.6" (1070 mm x 600 mm)
Height range	20.5" – 42.1" (520 mm – 1070 mm)	20.5" – 42.1" (520 mm – 1070 mm)	20.5" – 42.1" (520 mm – 1070 mm)
X-ray area	14.7" (370 mm)	Head end: 27.0" (685 mm) Foot end: 14.7" (370 mm)	Head end: 27.0" (685 mm) Foot end: 14.7" (370 mm)
Longitudinal shift	18.1" (460 mm)	15.7" (400 mm)	15.7" (400 mm)
Trendelenburg / reverse Trendelenburg	± 45°	± 45°	± 45°
Tilt	± 30°	± 30°	± 30°
Maximum patient weight	880 lbs (400 kg)	880 lbs (400 kg)	880 lbs (400 kg)

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.